

MR Paranasal Sinus Protocol

(** For Anosmia or CSF leaks please do "Anosmia/ CSF leak " protocol

| | | | | |
|---------------------------|-----|-----|-------|--|
| SAG T1 TSE | 4mm | | 22cm | |
| AX T2 FS | 3mm | 10% | 18-21 | Mouth opening to top of frontal sinuses |
| AX T1 | 3mm | 10% | 18-21 | Mouth opening to top of frontal sinuses |
| COR T2 | 3mm | 10% | 18-21 | Tip of nose to sella |
| COR T1 | 3mm | 10% | 18-21 | Tip of nose to sella |
| AX DW | 3mm | 10% | 18-21 | Same as axials (resolve/TSE DW if available) |
| <i>Contrast Injection</i> | | | | |
| AX T1 FS or VIBE/LAVA | 3mm | 10% | 18-21 | VIBE/ LAVA if available |
| COR T1 SE FAT SAT | 3mm | 10% | 18-21 | Same as pregad cor, no need for cor if VIBE/LAVA + |
| Ax BRAVO/MPGR BRAIN | 1mm | | | BRAIN, post only |
| | | | | |

