

Body MRI for Neurofibromatosis

Screening/Follow-up

(use small FOV if lesion specified; if not, use large FOV i.e. whole chest/body wall, extremity)

MRI Abd/Pelvis

Ax T2 (or COR T2)

Ax IP/OOP

Ax DWI

Ax/Sag/Cor STIR

Ax LAVA Pre

Ax/Cor/Sag LAVA Post

MRI Chest/Long Bone

Ax/Cor T1

Ax T1 FS

Ax T2

Ax/Sag/Cor STIR

Ax T1 FS Pre/Post

Cor/Sag T1 FS Post (or LAVA Post if not able to hold still)